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EN

Tampons are intended to be inserted into the vagina to absorb menstrual fluid. It is common for women using them for the first time to have to try several times to insert the tampon easily and comfortably.

1. Wash your hands and remove the wrapper from the tampon.
2. Assume a comfortable position, with your muscles relaxed. One of the positions that allows easy access to the vaginal area is standing up, with your knees slightly bent and separated, or with

one foot on the toilet.

3. Unravel the tampon string and pull it quickly and firmly. Make sure both ends of the string are connected.
4. Push your index finger into the hollow at the end of the tampon.
5. Insert the tampon into the vagina; this will be easier if the body is relaxed.
6. Push the tampon into the vagina as far back as possible, at an angle. The string will be hanging outside your body.
7. When the tampon is correctly placed, you will not feel anything. If you feel uncomfortable it means that you have not inserted the tampon to a sufficient depth. You should remove the tampon and insert a new one, ensuring that you insert it deeply.
8. Wash your hands.

WHEN TO CHANGE THE TAMPON

You should change your tampon every 4–6 hours. Wash your hands before and after removing the tampon. To remove the tampon, slightly bend your knees with your legs open and relax. Pull the string out in the same direction as you inserted it. Throw them away or dispose of them according to your country's legislation. If you feel discomfort when removing the tampon and can see white fibers, it may not be fully saturated yet. Always use an absorbency lower than your flow. Make sure you have always removed the last tampon before inserting another one. Ensure you remove the last tampon of the period and do not use tampons before or between menstruations.

ABSORBENCY CAPACITIES

6–9 g **Regular** for days with light to moderate flow
9–12g **Super** for days with moderate to heavy flow

WHAT YOU SHOULD KNOW ABOUT TOXIC SHOCK SYNDROME (TSS)

TSS is a rare but serious and potentially life-threatening illness. It is caused by toxins produced by *Staphylococcus aureus*, a bacterium commonly found on the skin, nose, and vagina. TSS can occur in women, men, and children. Menstrual TSS is linked to wearing an internal-use product for too long during menstruation. The risk of developing TSS is higher in teenagers and women under 30 years old than in older women. The estimated incidence of TSS is 1 to 17 cases per 100,000 women and girls of childbearing age. Although it is uncommon, it is important to be aware of it to know how to act if necessary.

Recognizing symptoms and early treatment are crucial. Symptoms that may indicate possible TSS include sudden high fever (102°F / 39°C or higher),

dizziness, vomiting, diarrhea, fainting or feeling faint when standing up, muscle aches, skin irritations resembling sunburn, or sore throat. Not all symptoms may be present, only some. In later stages of the illness, the skin may begin to peel. TSS can occur during menstruation or shortly after. It can progress rapidly with flu-like symptoms to a severe and potentially fatal illness. If you experience one or more of these symptoms, remove the tampon, keep it for analysis if possible, and consult your doctor immediately. Inform your doctor that you have been using tampons and that you are concerned about TSS. TSS may recur. If you have had TSS, the use of internal products such as tampons is not recommended. You should talk to your doctor before using tampons again or if you have questions about TSS or product use. You can reduce the risk of menstrual TSS by using tampons with the lowest absorbency needed to meet your needs and changing them regularly every 4–6 hours. You can also alternate an external product such as a sanitary towel or a menstrual pad instead of a tampon at least once a day. The use of tampons at night is discouraged; external products are recommended to reduce the risk of TSS. If you do not use tampons, the risk of menstrual TSS is reduced but not completely eliminated.

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